# INSTRUCTIONS FOR USE OF THE ONCOhabitats MEDICAL DEVICE









### HISTORIAL DE VERSIONES

Versión	Fecha
1.0	01/05/2023
2.0	12/07/2023

## INDEX

0 INTRODUCTION AND BACKGROUND INFORMATION	4
1 REGISTRATION IN THE SYSTEM	5
2 LOGGING ON TO THE SYSTEM	7
3 PERFORMING MORPHOLOGICAL SEGMENTATION	8
4 RUNNING A SEGMENTATION OF HABITATS VASCULAR	. 14
5 BIBLIOGRAPHY	22

### **0 INTRODUCTION AND BACKGROUND INFORMATION**

This document shows the instructions for using the ONCOhabitats software (patent number ES201431289A in Spain, EP3190542A1 in Europe, and US20170287133A1 in the USA). The product has been developed by the Biomedical Data Laboratory (BDSLab) of the Universitat Politécnica de València (UPV), Camino de Vera s/n 46022 València - Spain.

ONCOhabitats is a software for evaluating the heterogeneity of high-grade gliomas (IDHwild Glioblastoma and IDH-mutant Astrocytoma) using MRI data. The system offers two services: 1) segmentation of high-grade gliomas and 2) delineation of vascular heterogeneity of highgrade gliomas using the HTS method. From the input data of the case to be studied, provided by the user on the web page (https://www.oncohabitats.upv.es/), the system produces as output a report with the results obtained that is sent to the user by e-mail, in addition to displaying them on the web page. This information can be used to assist in understanding the physiology and anatomy of the tumor and thus to support physicians in planning surgeries within the framework of clinical research with medical devices. Current security standards will store all data (user-entered data and results generated by the software). As for end users, the software is open source but is intended for use by clinicians and researchers who wish to obtain and analyze information on high-grade gliomas.

The software was published on the website https://www.oncohabitats.upv.es/ for academic purposes on May 3, 2021. The current version of ONCOhabitats has been carried out within the ALBATROSS (PID2019-104978RB-I00) and SINUE (INNEST/2022/87) projects. It has been retrospectively evaluated in hospitals in Spain, Italy, Belgium, and Norway in Clinical Study NCT03439332 (Del Mar Álvarez-Torres et al., 2020). Currently, ONCOhabitats is being validated by the clinical study "Evaluation of the feasibility of ONCOhabitats to assist surgery and treatment planning in patients with Glioblastoma, IDH-wild type (INNEST/2022/87)". This clinical validation is being funded by the ALBATROSS project (PID2019-104978RB-100 - Agencia Estatal de Investigación) and the SINUÉ project (INNEST/2022/87 - Agencia Valenciana de la Innovación).

As a risk prevention action for a product under clinical validation, a training course will be conducted before recruiting the first patient, and multimedia information from the course will be available on the ONCOhabitats website to ensure adequate training of professionals on the use, interpretation, and limitations of the software concerning healthcare aspects.

Given the current innovative nature of AI-based technology, all AI-based software designed and developed for healthcare use should be considered by users as an aid in their analysis of patients and never to guide diagnoses or treatments on their own. However, it is globally acknowledged that specific risks are associated with data science and AI-based technology developments. Therefore, we must identify and assess the specific risks associated with our research projects, including potential vulnerabilities in data security and ethical or privacy issues.

The Biomedical Data Laboratory (BDSLab) of the Universitat Politécnica de València (UPV) should be notified in case of detecting any failure or malfunction of the product causing a severe incident. The contact methods are available in the "contact" section of the group's website: https://www.bdslab.upv.es/#contact-page. BDSLab, as a research group of the Universitat

Politècnica de València, promoter of the study, will notify the Agencia Española de Medicamentos y Productos Sanitarios (AEMPS) through the following e-mail: psinvclinic@aemps.es.

ONCOhabitats is hosted on a secure server at the UPV, under secure protocols and following the premises of the data protection delegate of the University itself. As for the necessary hardware, it is physically located in a server room with an electronic lock that can only be accessed by authorized personnel. Likewise, the data entered by the user and generated by the software are stored following the security standards in force.

## **1 REGISTRATION IN THE SYSTEM**

- 1.1. Accessing the website: <u>www.oncohabitats.upv.es</u>
- 1.2. Access the "REGISTER" section.



Figure 1. ONCOhabitats website.

1.3. Fill in the required data and click on the "captcha" as a test to check that the user is a person.

USERNAME	
usuario	
FIRST NAME	
nombre	
LAST NAME	
apellidos	
INSTITUTION	
institución	
E-MAIL ADDRESS	
nombre@dominio.com	
PASSWORD	
CONFIRM PASSWORD	
I'm not a robot	reCAPTCHA Privacy - Terms
Register	Login

Figure 2. Registration form.

After clicking on "Register" you will receive an email in the account indicated when registering, welcoming you. Your ONCOhabitats account will now be created.

## 2 LOGGING ON TO THE SYSTEM

- 2.1. Access the website: www.oncohabitats.upv.es
- 2.2. Access to the "LOGIN" section



Figure 3. ONCOhabitats website.

2.3. Enter your username or email account and password. Then click on "I'm not a robot" and, after confirming the CAPTCHA, click on the "Login" button.

USERNAME OR E-MAIL					
usuario					
PASSWORD					
•••••					
I'm not a robot	reCAPTCHA Privacy - Terms				
Login	Register				
Forgot your	password?				

Figure 4. Login screen.

## **3 PERFORMING MORPHOLOGICAL SEGMENTATION**

3.1. Accessing the website: <u>www.oncohabitats.upv.es</u>

3.2. Go to the section "SERVICES > SEGMENTATION SERVICE".



Figure 5. Selection of morphological segmentation service.

3.3. Access the "MR IMAGING" section. In the four boxes we introduce the MR image (T1, T1c, T2 and FLAIR) indicated in each one, either by dragging it from the folder where they are located, or by clicking on "BROWSE", looking for the folder where the image is located in the window that appears and, once located, click on Open.

O Abrir				×
← → ∽ ↑ 📙 « Pipe	linesShared > test > T1	~ Ō	,	
Organizar 👻 Nueva car	oeta			
<ul> <li>OneDrive</li> </ul>	Nombre		Fecha de modificación	Тіро
OneDrive - Persor	<b>X</b> T1		21/06/2022 16:51	Archivo NII
Este equipo Descargas Documentos Escritorio Elmágenes Música Objetos 3D Vídeos SC(C;) DATA (D:)				
				>
No <u>m</u> br	e: T1	~	Todos los archivos	~
			Abrir	Cancelar
T1 MRI BROWSE	het	<b>B</b> F	IC MRI	ł

Figure 6. Select MR image.

3.4. Fill in the fields in the "CLINICAL INFORMATION" section. The field "MRI storage consent:" indicates whether the user gives permission for the case entered into the platform to be stored on the platform. If permission is not given, the case in question will be deleted from the system some time after its execution. In a clinical trial, the images will not be saved on the platform's storage server, so this box should not be ticked.

## CLINICAL INFORMATION

RUN ANATOMIC	AL JOB
Edit	
MRI storage consent:	Help us to improve ONCOhabitats (Details) 🗹
Sex:	<ul> <li>Male</li> <li>Female</li> <li>Not Available</li> </ul>
Age:	50
ID:	

Figure 7. Clinical data entry.

3.5. Click on the "RUN ANATOMICAL JOB" button, which can be seen in the previous illustration. After this, you will see a screen indicating that the service has been successfully launched.



*Figure 8. Screen indicating the successful launch of the service.* 

After this, the system will start executing the case. Once finished, the system will send an email to the user with the results: on the one hand, a report in pdf format, and, on the other hand, a folder with the images generated during the execution.

ONCOhabitats job results 🕨 Recibidos 🗙
ONCOhabitats <info@oncohabitats.upv.es></info@oncohabitats.upv.es>
Dear ,
Your job test has been successfully processed!
The report and results can be downloaded from the following links:
Report Results
Thank you for using ONCOhabitats services,
ONCOhabitats
NOTE: For practical reasons, these results will only be accessible in our servers for the following 15 days
Figure 9. Mail received with the results.

You can also view the results by accessing the "MY JOBS" section of the website.



Figure 10. ONCOhabitats website.

Here you can view the data of the executed case and download the full report and results.

٢	2342	Anatomical segmentation	Finished	2023-01-02 11:58:38	2023-01-02 11:58:49	2023-01-02 12:37:33
		patient:			tes	t
gender:			other			
age:			50			
report:			download			
results:				down	load	

Figure 11. Window with the results of the job, shown in the "My Jobs" section.

The report generated and sent to the user follows the following structure:

-A cover page, with the study information data (report ID, report date, user name and institution and patient ID, gender and age), as shown in Figure 12.

A section with the calculated volumetry data for the morphological segmentation of highgrade gliomas (in cm<sup>3</sup>), the intracranial cavity mask (in axial, sagittal and coronal planes) and the tumor segmentation map (in axial, sagittal and coronal planes) generated during the execution of the service, as shown in Figure 12.

ONCOhabitats - Morphological Segmentation Report

	AMATION						
eport ID 2518e1f77	2166125ad3	d4c0aef3	37f7420e661ea		R. 0	eport Date 9/04/2021 - 10	:16:48
User Name jajuaal1			Institution Universitat Politecnica de Vale			E-mail jajuaal1@ibime.upv.es	
Patient ID AnatomicalSegmentationDemo			Patient Gender male			Patient Age 65	
BM MORPH	HOLOGICAL SI	EGMENT	ATION				
olumetry (	(cm <sup>3</sup> )						
Intra-cran	ial cavity	Enhand	cing Tumor	Edema		Necrosis	
Abs.	Rel.	Abs.	Rel.	Abs.	Rel.	Abs.	Rel.
1501.61	100.00%	25.55	1.70%	59.55	3.97%	13.81	0.92%
10000	Letter 1		11	2			
GBM anatol	mical segmen	ntation m	aps (axial / sag	ilttal / corona			

<sup>1</sup>The results contained in this report are only certified for research purposes.

Figure 12. Page 1 of the report generated by the morphological segmentation service.

- www.oncohabitats.upv.es -

A back cover showing disclaimer information, a definition of the processes carried out in the service (pre-processing of MR images and segmentation of gliomas of some grade) and a list of the acronyms used in the report, as seen in Figure 13.

ĩ

#### DISCLAIMER

All calculations, measurements and images provided by this software are intended only for scientific research. Any other use is entirely at the discretion and risk of the user. If you do use this software for scientific research please give appropriate credit in publications. The results of the HTS may not be commercially used in any other way without prior approval of the author.

#### DEFINITION OF THE SERVICES

**Pre-Processing:** Pre-processing module attempts to enhance and correct the MR images for its posterior analysis. Several common artifacts are corrected in this module such as magnetic bias field inhomogeneities, noise or motion artefacts. Additionally, automated registration and skull-stripping is conducted to generate a consistent high quality imaging data of the brain.

Segmentation: The anatomical analysis of the glioblastoma requires the delineation of the tumor tissues, which encompass the enhancing tumor, the edema and the necrotic tissues. In this module, we implemented a deep learning approach to provides such tissue identification. Convolutional Neural Networks are employed in combination with a morphological component analysis to identify the pathological structures.

- www.oncohabitats.upv.es -

ii -

*Figure 13. Back cover of the report generated by the morphological segmentation service.* 

## **4 RUNNING A SEGMENTATION OF HABITATS VASCULAR**

- 4.1. Accessing the website: <u>www.oncohabitats.upv.es</u>
- 4.2. Go to "SERVICES > VASCULAR HABITATS SERVICE".



Figure 14. Selection of vascular habitat segmentation service.

4.3. Access the "MR IMAGING" section. In the five boxes we introduce the MR image (T1, T1c, T2, FLAIR and DSC), indicated in each one either by dragging it from the folder where it is located, or by clicking on "BROWSE", looking for the folder where the image is located in the window that appears and, once located, click on Open.

O Abrir	×
$\leftarrow$ $\rightarrow$ $\checkmark$ $\uparrow$ $\bigcirc$ $\sim$ PipelinesShared $\rightarrow$ test $\rightarrow$ T1 $\checkmark$ $\bigcirc$	,> Buscar en T1
Organizar 🔻 Nueva carpeta	E 🕶 🔟 💡
OneDrive     Nombre	Fecha de modificación Tipo
T1	21/06/2022 16:51 Archivo NII
Este equipo Descargas Documentos Escritorio Escritorio Música Objetos 3D Vídeos G (C:) DATA (D:)	
No <u>m</u> bre: T1	Todos los archivos
	<u>A</u> brir Cancelar
T1 MRI BROWSE help help	F1C MRI

Figure 15. Select MR image

In addition, the user will have to enter the TE (Echo Time) and TR (Repetition Time) fields.

		 TE:	22,3	ms	
	DROP HERE <b>DSC MRI</b>	TR:	1500	ms	
	BROWSE	Edit			
help					



4.4. Fill in the fields in the "CLINICAL INFORMATION" section. The field "MRI storage consent:" indicates whether the user gives permission for the case entered into the platform to be stored on the platform. If permission is not given, the case in question will be deleted from the system some time after its execution. In the case of a clinical trial, the images will not be saved on the platform's storage server, so this box should not be ticked.

## CLINICAL INFORMATION

ID:	test
Age:	50
Sex:	<ul> <li>Male</li> <li>Female</li> <li>Not Available</li> </ul>
MRI storage consent:	Help us to improve ONCOhabitats (Details)
Edit RUN HTS JOB	
Figure	e 17. Clinical data entry.

4.5. Click on the "RUN HTS JOB" button, which can be seen in the previous illustration. After this, you will see a screen indicating that the service has been successfully launched.



Figure 18. Screen indicating the successful launch of the service.

After this, the system will start running the case. Once finished, the system will send an email to the user with the results attached: on the one hand, a report in pdf format, and, on the other hand, a folder with the images generated during the execution.

ONCOhabitats job results 🔉 Recibidos 🗙
ONCOhabitats <info@oncohabitats.upv.es></info@oncohabitats.upv.es>
Dear ,
Your job test has been successfully processed!
The report and results can be downloaded from the following links:
Report Results
Thank you for using ONCOhabitats services,
ONCOhabitats
NOTE: For practical reasons, these results will only be accessible in our servers for the following 15 days.
Figure 19. Mail received with the results.

You can also view the results by accessing the "MY JOBS" section of the website, as shown in Figure 20. Here you can see the data of the executed case and download the report and the complete results.

•	2341	Vascular habitats	Finished	2023-01-02 10:17:12	2023-01-02 10:17:49	2023-01-02 11:14:56		
		patient:		test				
gender:					othe	er		
		age:		50				
		TE:		null				
		TR:			nul	I		
		report:		download				
		results:			downl	oad		

*Figure 20. Window with the results of the job, shown in the "My Jobs" section.* 

The report generated and sent to the user follows the following structure:

-A cover page, with the study information data (report ID, report date, user name and institution and patient ID, gender and age), as shown in Figure 21.

A section with the calculated volumetry data for the morphological segmentation of highgrade gliomas (in cm<sup>3</sup>), the intracranial cavity mask (in axial, sagittal and coronal planes) and the tumor segmentation map (in axial, sagittal and coronal planes) generated during the execution of the service, as shown in Figure 21.

ONCOhabitats - Vascular Heterogeneity Assessment Report

5d9f1715e0	8472c3d79a	2ea298fdb	1ed7fb42c58	3	F	Report Date 09/04/2021	- 10:17	:49
User Name <i>jajuaal1</i>		li L	nstitution Universitat Po	olitecnica de	Vale j	-mail ajuaal1@ib	ime.up	1.es
Patient ID VascularHat	oitatsDemo	F /	Patient Gende male	er	F	Patient Age 55		
GBM MORPH	IOLOGICAL S	EGMENTAT	NOR					
Volumetry (	cm <sup>3</sup> )							
Intra-crani	al cavity	Enhancit	ng Tumor	Edema		Necr	rosis	
Abs.	Rel.	Abs.	Rel.	Abs.	Rel.	Abs.	8.	Rel.
1501.83	100.00%	25.47	1.70%	59.81	3.98%	6 14.0	14	0.94%
	17	))						5 /)
GBM anator	nical segmen	ntation map	Ds (axial / say	gittal / corona	al)			
GBM anator	nical segment	ntation map	os (axial / sag	giltal / corone				

*Figure 21. Page 1 of the report generated by the vascular habitats service.* 

-A section with the generated DSC quantification information, such as the median and mean absolute deviation of the DSC biomarkers per tissue (active tumor and edema), the biomarker maps (rCBV map, rCBF map and MTT map, in the axial plane) and the K2 leakage correction map, also in the axial plane, as seen in Figure 22.

#### ONCOhabitats - Vascular Heterogeneity Assessment Report

#### DSC QUANTIFICATION

- Median and MAD deviation of DSC biomarkers per tissue (enhancing tumor, edema and necrosis)

	rCBV	rCBF	MTT	
Enh.Tumor	3.07±0.88	2.24±0.66	1.29±0.11	
Edema	2.27±0.97	2.00±0.78	1.07±0.10	

- Blomarker maps: 1) rCBV map, 2) rCBF map, and 3) MTT map



- K2 leakage correction map (Boxerman et. al 2006)



results contained in this report are only certified for research purposes.	
- www.oncohabitats.upv.es -	н

Figure 22. Page 2 of the report generated by the vascular habitats service.

A section with the generated Hemodynamic Tissue Signature DSC perfusion maps and the metrics calculated in the service, such as volumetry per vascular habitat (in cm<sup>3</sup>), median and mean absolute deviation of DSC biomarkers per vascular habitat, the Hemodynamic DSC Tissue Signature nosology map (in axial, sagittal and coronal planes), the DSC Hemodynamic tissue signature chart and the prototype concentration-time perfusion curves, as shown in Figure 23.

The DSC Hemodynamic tissue signature chart shows how each DSC biomarker (CBV, CBF, MTT and K2) calculated in the DSC quantification process relates to each habitat (HAT, LAT, IPE and VPE), and the prototype concentration-time perfusion curves show the evolution of the intensity of the bolus of contrast injected into the patient in each habitat over time.

HEMODYN/	MIC DSC TIS	SUE SIGNAT	URE				
- Volumetry	per vascular	habitat (cm	3)				
HAT				IPE		VPE	
Abs.	Rel.	Abs.	Rel.	Abs.	Rel.	Abs.	Rel.
11.94	0.80%	25.52	1.70%	11.05	0.74%	20.26	1.35%

- Median and MAD deviation of DSC biomarkers per vascular habitat

25.52

	rCBV	rCBF	MTT	
HAT	5.16±0.61	4.30±0.67	1.21±0.10	
LAT	3.25±0.45	2.47±0.50	1.32±0.13	
IPE	2.07±0.30	1.72±0.23	1.21±0.07	
VPE	0.97±0.28	0.99±0.20	0.96±0.10	

- Hemodynamic DSC tissue signature nosological map



- Hemodynamic DSC tissue signature chart and prototype concentration-time perfusion curves



*Figure 23. Page 3 of the report generated by the vascular habitats service.* 

A back cover showing disclaimer information, a definition of the processes carried out in the service (pre-processing of MR images and segmentation of gliomas of some grade) and a list of the acronyms used in the report, as seen in Figure 24.

#### DISCLAIMER

All calculations, measurements and images provided by this software are intended only for scientific research. Any other use is entirely at the discretion and risk of the user. If you do use this software for scientific research please give appropriate credit in publications. The results of the HTS may not be commercially used in any other way without prior approval of the author.

#### DEFINITION OF THE SERVICES

**Pre-Processing:** Pre-processing module attempts to enhance and correct the MR images for its posterior analysis. Several common artefacts are corrected in this module such as magnetic bias field inhomogeneities, noise or motion artefacts. Additionally, automated registration and skull-stripping is conducted to generate a consistent high quality imaging data of the brain.

Segmentation: The anatomical analysis of the glioblastoma requires the delineation of the tumor tissues, which encompass the enhancing tumor, the edema and the necrotic tissues. In this module, we implemented a deep learning approach to provides such tissue identification. Convolutional Neural Networks are employed in combination with a morphological component analysis to identify the pathological structures.

DSC perfusion quantification: DSC quantification involves the computation of the hemodynamic indices obtained from a kinetic analysis of the first pass of a intravenously injected paramagnetic contrast agent. T1 and T2 leakge effects are also corrected to not miss-estimate the hemodynamic biomarkers. The quantified maps computed by this module are: relative Cerebral Blood Volume (rCBV), relative Cerebral Blood Flow (rCBF), Mean Transit Time (MTT) and K2 maps.

Hemodynamic Tissue Signature (HTS): HTS consist on an automated unsupervised method able to describe the vascular heterogeneity of the enhancing tumor and edema tissues in terms of the angiogenic process located at these regions. The HTS provides a characterization of the GBM, whose output is a nosologic map of the tumoral tissues grouped in different vascular sub-compartments with their associated MRI fingerprint. We consider 4 vascular sub-compartments for the GBM, the high angiogenic enhancing tumor region (HAT), the low angiogenic enhancing tumor region (LAT), the potentially tumor infiliatrated peripheral edema (IPE) and the pure vasogenic edema (VPE). The HTS is able to capture the local heterogeneity of the tumor, hence providing relevant information about its behaviour.

rCBV relative Cerebral Blood Volume

- rCBF relative Cerebral Blood Flow
- MTT Mean Transit Time
- DSC Dynamic Susceptibility Contrast
- HTS Hemodynamic Tissue Signature

- www.oncohabitats.upv.es -

iv

Figure 24. Back cover of the report generated by the vascular habitats service.

## **5 BIBLIOGRAPHY**

Del Mar Álvarez-Torres, M. *et al.* (2020) 'Robust association between vascular habitats and patient prognosis in glioblastoma: An international multicenter study', *Journal of magnetic resonance imaging: JMRI*, 51(5), pp. 1478–1486. Available at: https://doi.org/10.1002/jmri.26958.